Quartet Data Portal User Guide for

Proteomic Quality Assessment Report


1. Preparation

QDP Address	https://chinese-quartet.org/
Input files	The file format (suffix) should be .csv.

2. Log on to QDP

2.1 QDP account


1) If you have a request for the reference materials, please go to http://chinese-quartet.org/#/materials to request the reference materials. The Quartet team will contact you for further confirming your information and send you a registration email.




2) If you do not have a request for the reference materials, please send an email to **quartet@fudan.edu.cn** for a QDP account.


2.2 Log on to QDP

Address: http://chinese-quartet.org/






2.3 Select Quality Assessment> QC Apps





3. Run QC Report


3.1 Select "Proteomics" > "Run".



3.2 Click the symbol called "New QC Report".









3.4 Select "Step2: Upload File(s)" > "QC Report for Quartet Proteomics". Please note that the files need to be in csv format. For details of the fixed format of the uploaded files, please refer to

https://docs.chinese-quartet.org/data_pipelines/proteomics/qc_report/.



3.5 After uploading the tested files, select "Step3: Parameters & Submit". Fill in the required information and click "Submit".



3.6 Then you can preview or download your QC report by clicking "Result" or "Download" respectively.

